


ANALYTICAL CONDITIONS

INSTRUMENT	SHIMADZU LC-20AT HPLC SHIMADZU SPD-20A UV DETECTION		
COLUMN	GL Sciences InertSustain C18 4.6 x 150 mm, 5µm particle size; S/N 3AR97176		
MOBILE PHASE	A – Methanol, B –0.1mol/L Ammonium acetate; A:B=65:35	Product Name:	Cyclopamine
COLUMN TEMP.	30 °C	InvivoChem Cat #:	V1332
FLOW RATE	1.0 ml/min	CAS #:	4449-51-8
INJECTION VOL.	5 µl	Batch #:	V133202
INJECTION CONC.	2 mg/ml in MeOH	Formula:	C27H41NO2
DETECTION	UV 210 nm	Molecular Weight:	411.62



分析结果表

峰号	峰名	保留时间	峰高	峰面积	含量
1		6.388	4861.918	100668.203	1.2682
2		9.918	225708.203	7836977.500	98.7318
总计			230570.122	7937645.703	100.0000