

File D:\DATA

Injection Date : 1 Feb 16 2:45

Sample Name : 841290-80-0

Acq. Method : D:\METHODS\1-POS-3MIN.M

Analysis Method : D:\METHODS\1-POS-3MIN.M

Sample Info : Easy-Access Method: '1-POS-3min'

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)

Gradient: 5%-95% B in 1.2min


Flow Rate: 2.2 ml/min

Column: Poroshell 120 EC-C18, 4.6*30mm, 2.7um


Oven Temperature: 50 C

Mass Range: 110-1500


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.44	0.02	61.78	46.87	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.39	0.02	17.37	16.67	0.50
1.44	0.02	3434.45	2896.29	98.92
1.70	0.01	20.20	24.55	0.58


Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.39	0.01	10.66	11.47	0.38
1.44	0.02	2776.79	2811.98	99.10
1.66	0.02	5.17	5.28	0.18
1.70	0.01	9.34	11.38	0.33


LC/MS Report

Ret. Time: 0.23


<<<< POSITIVE SPECTRA >>>>



Ret. Time: 0.38



Ret. Time: 1.45



Ret. Time: 1.71

