

HPLC Analysis for V5119401 RMC-6236


=====

Acq. Operator : SYSTEM Seq. Line : 14
 Acq. Instrument : 1260-DAD Location : P1-D10
 Injection Date : 6/17/2024 3:50:47 PM Inj : 1
 Inj Volume : 20.000 µl

Different Inj Volume from Sample Entry! Actual Inj Volume : 5.000 µl

Acq. Method : D:\DATA\2024\05\20240613001 2024-06-17 11-44-35\12MIN5-95.M
 Last changed : 6/17/2024 11:44:35 AM by SYSTEM
 Analysis Method : D:\DATA\2024\05\20240613001 2024-06-17 11-44-35\12MIN5-95.M (Sequence Method)
 Last changed : 6/17/2024 5:36:48 PM by SYSTEM
 (modified after loading)

Method Info : Column: Welch Ultimate XB- C18 4.6X150mm 5µm
 Mobile phase: A:10mmol NH4HCO3 B: ACN Flow: 1.0mL/min
 Time: 0.0 3.0 8.0 12.0
 B%: 5% 5% 95% 95%



=====

Area Percent Report

=====

Sorted By : Signal
 Multiplier : 1.0000
 Dilution : 1.0000
 Do not use Multiplier & Dilution Factor with ISTDs

Product Name: RMC-6236
 InvivoChem Cat #: V51194
 CAS #: 2765081-21-6
 Batch #: V5119401
 Formula: C44H58N8O5S
 Molecular Weight: 811.06

Signal 1: DAD1 B, Sig=254,4 Ref=off

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	7.150	BB	0.0501	5.03597	1.54070	0.0323
2	7.778	BBA	0.1640	16.54148	1.58191	0.1062
3	9.089	BBA	0.0720	15.88860	3.31055	0.1020
4	9.353	BBA	0.0631	18.69392	4.44022	0.1201
5	10.894	BBA	0.2558	1.55140e4	879.53308	99.6393

Totals : 1.55702e4 890.40647


*** End of Report ***

HPLC Analysis for V5119401 RMC-6236

SAMPLE INFORMATION

Sample Name: RMC6236
Sample Type: Unknown
Vial: 1:B,3
Injection #: 1
Injection Volume: 1.00 µl
Run Time: 15.0 Minutes
Date Acquired: 9/12/2024 09:11:40 AM CST
Date Processed: 9/12/2024 10:31:06 AM CST

Acquired By: ZHB
Sample Set Name: 20240912_116_01
Acq. Method Set: Hex_EtOH_60_40
Processing Method: RMC6236
Channel Name: 254.0nm
Proc. Chnl. Descr.: 2998 PDA 254.0 nm



	RT	Area	% Area	Height
1	3.510	170424	0.62	7465
2	5.647	27160701	99.38	690853

Product Name: RMC-6236
InvivoChem Cat #: V51194
CAS #: 2765081-21-6
Batch #: V5119401
Formula: C44H58N8O5S
Molecular Weight: 811.06