

HPLC Analysis for V423201 Hibifolin




Chengdu Gelipu Biotech Co Ltd

SAMPLE INFORMATION

Sample Name: 棉黄素-8-0PTTG 07
Sample Type: Standard
Vial: 16
Injection #: 1
Injection Volume: 10.00 ul
Run Time: 100.0 Minutes
Date Acquired: 2024-5-10 9:54:12 CST
Date Processed: 2024-5-10 10:49:23 CST

Acquired By: System
Sample Set Name:
Acq. Method Set: 棉黄素80PTTG
Processing Method: Quality Guarantee System2
Channel Name: 254.0nm
Proc. Chnl. Descr.: PDA 254.0 nm



| | RT | Area | % Area | Height |
|---|--------|---------|--------|--------|
| 1 | 4.344 | 19784 | 0.21 | 1638 |
| 2 | 4.930 | 26209 | 0.28 | 3027 |
| 3 | 6.167 | 33945 | 0.36 | 3180 |
| 4 | 10.025 | 9296280 | 99.02 | 676635 |
| 5 | 23.120 | 12142 | 0.13 | 2255 |

Reported by User: System
Report Method: Chengdu Gelipu Biotech Co
Report Method ID: 124299
Page: 1 of 1

Project Name: Defaults
Date Printed: 2024-5-10
10:49:50 PRC