



Sample Information

Product Name : TY-A-10034 Pinealon
 Catalog ID : TY-A-10034 Pinealon
 Lot. No : P26031510034
 Sequence : Pinealon (H-Glu-Asp-Arg-OH)
 Column : Agilent 5 HC-C18(2) 250*4.6mm
 Buffer : A: 0.1% TFA in Acetonitrile B: 0.1% TFA in H2O
 Gradient : 2-27% in25min
 Flow rate : 1.0ml/min
 Wavelength : 220nm
 Dissolution method : 100% H2O
 Column temperature : 35°C



RT	Area	% Area	Height
5.643	9392066	99.11	727356
6.223	53883	0.57	8083
7.127	30149	0.32	6067



Product Name: Pinealon
 InvivoChem Cat #: V69561
 CAS #: 175175-23-2
 Batch #: V6956101
 Formula: C15H26N6O8
 Molecular Weight: 418.40