


## High Performance Liquid Chromatography

NO:023-K81-TSX

<b>Product Name</b>	Kisspeptin-10
<b>Sequence</b>	Tyr-Asn-Trp-Asn-Ser-Phe-Gly-Leu-Arg-Phe-NH <sub>2</sub>
<b>Dissolution method</b>	100%H <sub>2</sub> O
<b>Gradient</b>	15-45% in 25 min
<b>Flow rate</b>	1.0ml/min
<b>Wavelength</b>	220nm
<b>Column</b>	COSMOSIL Packed Column,5C18-MS-II,(4.6ID . 250mm)
<b>Buffer</b>	A: 0.1% TFA in Acetonitrile    B: 0.1% TFA in H <sub>2</sub> O



	RT	Area	%Area	Height
1	17.688	18090	0.34	1719
2	18.158	5190903	98.72	609335