

HPLC Analysis for ML-385



Method Info : Column:XBridge C18 (4.6x 50 mm, 3.5um)

Mobile phase: H2O(10 mmol NH4HCO3) (A) / ACN(B)

Elution program: Gradient from 10 to 95% of B in 1.5min at 1.8ml/min

Temperature: 50 degrees

Detection: UV (214 , 4 nm) and MS (ESI, POS mode)

Product Name: ML-385

InvivoChem Cat #: V2828




CAS #: 846557-71-9

Batch #: V282802

Formula: C29H25N3O4S

Exact Mass: 511.16

Molecular Weight: 511.60



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.82	0.01	839.03	835.11	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.82	0.01	337.63	346.86	100.00

Ret. Time: 0.56

<<<< POSITIVE SPECTRA >>>>

*ES-API Positive <<< Below Cutoff >>>

Ret. Time: 1.81



Ret. Time: 2.64

