



File D:\DATA

Injection Date : 12-Mar-19, 12:21:23

Sample Name : 700874-72-2

Acq. Operator : HPLC

Spec. Reported : MS Integration

Acq. Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M

Analysis Method : C:\Chem32\1\METHODS\2-HPLC-15MIN.M

Sample Info :

Method Info : Mobile Phase: A: water(0.01%TFA) B: ACN(0.01%TFA)

Gradient: 0min 5%B,3min 5%B,10min 95%B,15min 95%B

Flow Rate: 1.2 ml/min

Column: Eclipse XDB-C18,4.6*150mm,5um

Oven Temperature: 40 C


Tgt Mass(EZX) :

Location : P1-B-04


Inj : 0

Inj Volume : 1 ul

DAD1 A, Sig=214,4 Ref=off



DAD1 B, Sig=254,4 Ref=off



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
6.13	0.03	6.70	3.64	0.08
6.65	0.03	8140.60	3373.93	99.64
6.81	0.01	6.00	7.56	0.07
7.56	0.03	7.81	3.92	0.10
7.68	0.03	8.50	4.06	0.10

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
6.66	0.03	6617.94	3000.92	99.89
6.81	0.02	7.12	6.92	0.11