


HPLC Analysis for V006101 ISRIB (trans-isomer)


SAMPLE INFORMATION

Sample Name: ISRIB
Sample Type:
Acq. Method Set: Z
Injection #: 11
Injection Volume: 20.00 ul
Run Time: 30.0 Minutes
Proc. Chnl. Descr: 283

Acquire By: System
Date Acquired 2020-04-25 16:26:40
Date Processed: 2020-04-25 16:57:35
Processing Method: 283



	RT	Area	Height	% Area
1	16.374	43503	4670	0.10
2	17.340	41285382	3015941	99.07
3	17.989	97965	8422	0.23
4	18.676	311408	22128	0.24
5	19.645	326421	28261	0.28
6	21.149	34758	3964	0.08



Product Name: ISRIB (trans-isomer)
InvivoChem Cat #: V0061
CAS #: 1597403-47-8 (trans-ISRIB)
Batch #: V006101
Formula: C₂₂H₂₄Cl₂N₂O₄
Molecular Weight: 451.34