


(modified after loading)

DAD1 A, Sig=210,16 Ref=360,100 (110811\11081102.D)



=====
Area Percent Report
=====

Sorted By : Signal
Multiplier : 1.0000
Dilution : 1.0000

Signal 1: DAD1 A, Sig=210,16 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	10.990	BB	0.4581	3.57054e4	1156.91882	99.9483
2	13.190	BB	0.1263	7.30371	8.01367e-1	0.0204
3	14.818	BB	0.1507	11.18079	1.11804	0.0313

Totals : 3.57239e4 1158.83823

Results obtained with enhanced integrator!

=====
*** End of Report ***