




HPLC Analysis for V956201 4-Azido-L-phenylalanine

LC/MS Report

Acq. Method : D:\METHODS\1-POS-3MIN.M
Analysis Method : D:\METHODS\1-POS-3MIN.M
Sample Info : Easy-Access Method: '1-POS-3MIN'
Method Info : Mobile Phase : A:water(0.01%TFA) B:ACN(0.01%TFA)
Gradient:5%B for 0.2min,increase to 95%B within 1.5min,95%B for 1.5min,back within 0.01min
Flow Rate :2ml/min
Column :Sunfire, 50*4.6mm,3.5um
Column Temperature:50C



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
1.39	0.02	436.20	360.24	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.39	0.02	616.81	599.81	100.00



Product Name: 4-Azido-L-phenylalanine
InvivoChem Cat #: V9562
CAS #: 33173-53-4
Batch #: V956201
Formula: C9H10N4O2
Molecular Weight: 206.205

LCMS Analysis for V956201 4-Azido-L-phenylalanine

LC/MS Report

Ret. Time: 1.40

<<<< POSITIVE SPECTRA >>>>



Product Name: 4-Azido-L-phenylalanine
InvivoChem Cat #: V9562
CAS #: 33173-53-4
Batch #: V956201
Formula: C₉H₁₀N₄O₂
Molecular Weight: 206.205