


LCMS Analysis for V3639801 Cecropin P1

MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2023-12-16
Analyst: SHAO
Sample: 601131
M.W.: 3338.9
Lot. No.:

Instrument

Waters 2695
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Probe Bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.5ml/min
B. Conc.: 50%H2O/50%ACN