


MS Analysis for V662701 Examorelin

MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2024-12-11
Analyst: PY
Sample: DZT-020
M.W.: 887.04
Lot. No.: 20240922-2

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN