


MS Analysis for V3334301
C-Type Natriuretic Peptide (CNP) (1-22), human

MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2024-05-09
Analyst: SHAO
Sample: 601402
M.W.: 2197.6
Lot. No.:

Instrument

Waters 2695
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Probe Bias:

+4.5kv
Detector: 1.5kv
T. Flow: 0.5ml/min
B. Conc.: 50%H₂O/50%ACN