


LCMS Analysis for V6000501 Follicle-stimulating hormone

MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2023-01-05

Analyst: SHAO

Sample: 605130

M.W.: 980.2

Lot. No.:

Instrument

Probe:

Nebulizer Gas Flow:

CDL:

CDL Temp.:

Block Temp.:

Waters 2695

ESI

1.5L/min

-20.0v

250 °C

200 °C

Probe Bias:

+4.5kv

Detector:

1.5kv

T. Flow:

0.5ml/min

B. Conc.:

50%H2O/50%ACN