

1.31/1
243984-11-4

NMR Analysis for V459801 Resatorvid: Solvent: DMSO

^1H NMR (400 MHz, DMSO- d_6) δ 9.76 (s, 1H), 7.54 (dt, $J = 8.5, 3.9$ Hz, 2H), 7.27 (td, $J = 8.5, 2.9$ Hz, 1H), 7.10 (t, $J = 3.8$ Hz, 1H), 4.30 (d, $J = 5.6$ Hz, 1H), 4.00 (q, $J = 7.1$ Hz, 2H), 2.46 (d, $J = 13.7$ Hz, 1H), 2.33 (t, $J = 5.2$ Hz, 1H), 2.22 (s, 1H), 2.09 (d, $J = 9.2$ Hz, 1H), 1.76 (td, $J = 13.7, 12.8, 6.8$ Hz, 1H), 1.63 (d, $J = 11.7$ Hz, 1H), 1.05 (t, $J = 7.1$ Hz, 3H).

