

3.27/4
53-84-9 01batch

NMR Analysis for V3023303 NAD⁺; Solvent: D2O

¹H NMR (400 MHz, Deuterium Oxide) δ 9.33 (s, 1H), 9.18 (d, *J* = 6.3 Hz, 1H), 8.84 (dt, *J* = 8.2, 1.4 Hz, 1H), 8.51 (s, 1H), 8.31 (s, 1H), 8.19 (dd, *J* = 8.1, 6.3 Hz, 1H), 6.07 (dd, *J* = 10.3, 5.4 Hz, 2H), 4.56 – 4.40 (m, 4H), 4.38 (dd, *J* = 5.1, 2.6 Hz, 1H), 4.30 (h, *J* = 3.2 Hz, 2H), 4.20 – 4.08 (m, 3H).

