

NMR Analysis for V081001 Diosmetin; Solvent: DMSO

¹H NMR (400 MHz, DMSO-d₆) δ 12.94 (s, 1H), 10.86 (s, 1H), 9.46 (s, 1H), 7.55 (dd, *J* = 8.5, 2.3 Hz, 1H), 7.43 (d, *J* = 2.3 Hz, 1H), 7.09 (d, *J* = 8.6 Hz, 1H), 6.76 (s, 1H), 6.47 (d, *J* = 2.1 Hz, 1H), 6.20 (d, *J* = 2.1 Hz, 1H), 3.87 (s, 3H).

