


9708



```

NAME          Sep27-2013
EXPNO         50
PROCNO        1
Date_         20130927
Time_         14.57
INSTRUM       spect
PROBHD        5 mm PABBO HB-
PULPROG       zg30
TD            65536
SOLVENT       DMSO
NS            2
DS            2
SWH           8276.146 Hz
FIDRES        0.126314 Hz
AQ           3.9584243 sec
RG            181
DM           60.400 usec
DE           4.50 usec
TE           297.2 K
D1           1.00000000 sec
TD0          1
  
```

```

----- CHANNEL f1 -----
NUC1          1H
P1            12.58 usec
PL1           0.00 dB
PL1M          10.87646866 M
SFO1          400.1324710 MHz
SI            32768
SF            400.1300000 MHz
WDW           EM
SSB           0
LB            0.30 Hz
GB            0
PC            1.00
  
```