



8.199  
7.679  
7.659  
7.488  
7.468  
7.396  
7.393  
7.375  
7.358  
7.354  
7.324  
7.322  
7.303  
7.286

4.367

3.689  
3.678

3.390  
3.374  
3.330

2.740  
2.513  
2.508  
2.504  
2.499  
2.495

1.720  
1.700  
1.682  
1.662

1.346  
1.329  
1.311  
0.870  
0.852  
0.833

-0.000

1.91  
1.00  
0.99  
1.06  
1.08

2.08

2.00

3.72

2.10

2.16

8.86

3.14

9.0 8.5 8.0 7.5 7.0 6.5 6.0 5.5 5.0 4.5 4.0 3.5 3.0 2.5 2.0 1.5 1.0 0.5 0.0 -0.5  
f1 (ppm)