

2023-4-9/4
110505-75-4

^1H NMR (400 MHz, $\text{DMSO}-d_6$) δ 11.99 (s, 1H), 3.62 (dq, $J=8.8, 6.2$ Hz, 1H), 3.25 – 3.06 (m, 2H), 2.42 (dq, $J=12.4, 6.3$ Hz, 1H), 2.21 (t, $J=7.3$ Hz, 2H), 1.87 (dq, $J=13.4, 6.8$ Hz, 1H), 1.68 (dtd, $J=13.4, 8.3, 5.7$ Hz, 1H), 1.61 – 1.47 (m, 3H), 1.45 – 1.31 (m, 2H).

